

# COMPACT DRAINAGE PUMP

## Model 270



**MODEL 270**

**DRAINAGE UNIT** – The Gomco® Model 270 is designed for gastrointestinal drainage in the home, nursing homes, clinics or hospitals.

### FEATURES

- Lightweight for easy portability
- Ideal for homecare use
- Enclosed cabinet protects internal parts
- Low maintenance
- Hospital grade electric plug
- Two levels of drainage (90 mm Hg & 120 mm Hg)
- Positive shut-off collection system
- Bacteria filter
- Illuminated on/off switch and pump cycle light
- Floor bottle basket

**MODEL 270**– This Gomco unit is engineered to provide suction by alternating the expansion and contraction of air within a cylinder at regular intervals. There are no moving parts, resulting in a quieter, more reliable operation. The unit is designed to operate for an extended period of time, while solid state controls help to assure mild suction. This drainage unit includes two preset drainage levels (90 mm Hg and 120 mm Hg).

The bacteria filter reduces environmental hazards by not allowing pathogens to enter exhaust gas, and helps to prevent fluid and aerosol contamination of the unit. Unit includes a 1/2-gallon glass or disposable collection bottle, floor bottle carrier and two 6-foot lengths of disposable tubing. Model 270 is UL listed.

MODEL	COLLECTION BOTTLE	COMPRESSOR	VACUUM SETTING	FLOW RATE	DIMENSIONS HxDxW	SHIPPING WEIGHT
270	Glass	Thermotic Pump	90 mm Hg and 120 mm Hg	.25 lpm	11 x 8 x 12 in.	15 lbs.
2701	Disposable			.25 lpm	28 cm x 20 cm x 30 cm	7 kg

### ORDERING INFORMATION

w/Glass	w/Disposable	Voltage
01-21-0270	01-21-2701	110/50
01-22-0270	01-22-2701	115/60
01-73-0270	01-73-2701	220/50-60