



Vacuum traps help safeguard vacuum regulators, wall outlets and the central vacuum system from accidental contamination from canister overflows. This is accomplished by a mechanical float which cuts off suction in the event of a canister overflow.

**Key Features**

- Multiple inlet connectors are available to connect the vacuum trap to any vacuum regulator
- Standard hose barb outlet
- Reusable design can be used over and over

**Ordering Information**

<b>12-20-1000</b>	Vacuum Trap with DISS nut & nipple inlet
<b>12-20-2000</b>	Vacuum Trap with vacuum trap male fitting inlet
<b>12-20-3000</b>	Vacuum Trap with DISS hand tight inlet