




# Introducing the Schuco® Twist-Lock HD Collection Canister

The Schuco Twist-Lock HD single patient use collection canister features a convenient twist-off lid that can be effortlessly yet securely attached and removed without the potential for leaking or cracking that can occur with a snap-on lid. The dishwasher safe Twist-Lock canister is made of heavy duty 3 mm polycarbonate; the polypropylene lid has a silicone gasket seal. Clearly labeled patient and vacuum ports make assembling the canister a breeze. Schuco Twist-Lock HD canisters and lids are individually bagged and available in cases of 5 and 70.



**ORDERING INFORMATION**  
**S1160BS-RPL** Canister, 800 cc, 5/case  
**S1160BS-CS** Canister, 800 cc, 70/case



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
 314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.  
 Toll Free: 800/444-3940  
 FAX: 800/477-7701

Canada & Mexico  
 Toll Free: 800/446-0552  
 FAX: 800/246-6201

International  
 Telephone: 314/268-1683  
 FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 23-00-0007  
 Rev. June 2011