

LIFE+LITE™ Disposable Emergency Stretcher



The LIFE+LITE Disposable Emergency Stretcher is used for patient transport in virtually any situation. Constructed of strong corrugated plastic and weighing only 4.5 pounds, the LIFE+LITE can support up to 500 pounds.

Because of its construction and low cost, the LIFE+LITE can be easily disposed of, recycled, or cleaned and used again.

- Weighs only 4.5 lbs. (2 kg)
- Supports a patient up to 500 lbs. (227 kg)
- 10 large hand-holds for easy carrying
- Flexible PVC and corrugated plastic construction (not intended for spinal immobilization)
- Ideal for disaster management and preparation, as well as HazMat and patient transport for hospital, prehospital and industrial applications

SPECIFICATIONS

Length: 72" (183 cm)
Width: 22" (56 cm)
Depth: 1" at handle (2.5 cm)
Weight: 4.5 lbs. (2 kg)
Capacity: 500 lbs. (227 kg)
Material: Corrugated plastic/PVC

ORDERING INFORMATION

L724005 LIFE+LITE Emergency Stretcher
Sold only in cases of 5



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and
subject to change without notice.

Form 90-00-0022
June 2013