

# *Schuco S5000 Nebulizer*



*Introducing the Schuco® S5000 Nebulizer ... Combining Schuco quality and reliability with an affordable price.*

- Suppressed sound
- Enhanced airflow
- Easy to use and clean
- Our most cost-effective nebulizer

*Complete with all that's needed to begin nebulizer therapy:*

- Compact compressor
- Nebulizer cup
- Angled mouthpiece
- Oxygen tubing
- Five extra compressor filters
- 5-year limited warranty

**Schuco**  
byAllied

# Schuco<sup>®</sup> S5000

## STANDARD FEATURES

- Quiet operation ... less than 60 dBA
- Compact and lightweight ... just 3-1/2 lbs.
- Easy to use and clean
- Convenient cord storage compartment
- Substantial compressor pressure ... 35-45 psig
- High compressor operating flow rate ... 7-11 lpm
- Therapeutic particle size ... 0.5 µm-5µm
- AC operation ... 120V 60 Hz
- 5-year limited warranty



61399



61400

## ORDERING INFORMATION

**S5000** Compressor, nebulizer cup, angled mouthpiece, oxygen tubing, 5 extra compressor filters

## ACCESSORIES

**S66-0126** Replacement Filters, 5/pkg. (20 pkgs./case)

**64085** Adult Aerosol Mask with nebulizer and 7' clear tubing (50/case)

**64095** Pediatric Aerosol Mask with nebulizer and 7' clear tubing (50/case)

**61399** Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee, 6" flex tube (50/case)

**61400** Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee (50/case)

**88100** AC to DC Power Converter for use in car



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA  
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 23-00-0005  
Rev. Jan 2018