

*OptiVac<sup>®</sup>*  
*AC/DC Portable Suction Unit*  
*The Quality of Gomco<sup>®</sup> Where You Need It Most.*



The OptiVac is the most versatile product in the industry-leading Gomco line of portable suction equipment. Its compact design, ease of operation, rechargeable battery, AC power supply, and fully-enclosed unit are just a few of its unique features.

**Standard Features**

- 3 hours run time at maximum vacuum on a fully-charged battery
- 1 hour run time after just 2 hours of charging
- AC or DC operation ... ready when you are, 24/7
- 25-550 mm Hg vacuum range, coupled with greater than 30 lpm flow, provides powerful suction for faster aspiration
- Fully chargeable in 8 hours
- Will run at full performance on AC power even with a depleted battery
- Battery maintenance system keeps the battery fully charged but not over-charged
- Long-life internal rechargeable battery eliminates the need for frequent replacement
- Sealed water-resistant unit is easy to clean, and minimizes potential contamination or damage to the compressor and battery.
- Angled control panel for easy viewing
- Impact-resistant and durable
- 3-year warranty on the unit, 1 year on the battery

# OptiVac® AC/DC Portable Suction Unit



Tubing can be cut to desired length and attached to disposable hydrophobic/bacteria filter.



Carrying Case



Disposable Collection Canister



Disposable Hydrophobic/  
Bacteria Filter



Suction Tubing Kit

## SPECIFICATIONS

<b>Dimensions:</b>	9.4" H x 7.5" W x 16.8" L (239 mm x 191 mm x 427 mm)
<b>Weight:</b>	11.4 lbs. (5.18 kg)
<b>Vacuum Range:</b>	25-550 mm Hg
<b>Vacuum Flow:</b>	Greater than 30 lpm
<b>AC Operation:</b>	110-230 volts, 50/60 Hz
<b>DC Operation:</b>	12 volt DC battery
<b>Battery Type:</b>	Sealed lead acid
<b>Battery Life at Maximum Vacuum:</b>	3 hours at 550 mm Hg
<b>Battery Charge Time:</b>	Full charge in 8 hours
<b>Battery Charge Status:</b>	Yes; 5 indicator lights
<b>Bad Battery Sensor:</b>	Yes; 5 red LEDs
<b>Battery Depleted:</b>	Unit automatically shuts off
<b>AC Charger:</b>	Internal charger
<b>Operation with Depleted Battery:</b>	Yes; charging circuit can run on AC if battery is depleted
<b>Operating Temperature:</b>	-4° F to 120° F
<b>Sound Level:</b>	60 dBA
<b>UL Approval:</b>	Tested to UL 60601-1
<b>Warranty:</b>	3 years on unit, 1 year on battery

## ORDERING INFORMATION

<b>G180</b>	OptiVac® Portable Suction Unit with gauge & regulator, ON/OFF buttons, 1500 ml disposable canister, suction tubing kit, 3 bacteria filters, 115/60
<b>G180CE</b>	OptiVac® Portable Suction Unit; 230/50

## ACCESSORIES

<b>L190-BAG</b>	Carrying Case with shoulder strap
<b>20-08-0003</b>	Disposable Collection Canister, 1500 ml, stem inlet (48/case)
<b>20-08-0016</b>	Disposable Collection Canister, 1500 ml, stem inlet (16/case)
<b>01-90-3928</b>	Disposable Hydrophobic/Bacteria Filter, 1/4" hose barb x 1/4" hose barb (3/pkg.)
<b>S610100</b>	Suction Tubing Kit: 1 each 72" & 15" length, 1/4" ID, blue-tipped, clear PVC



**ISO 13485: 2016**

Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110  
314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • www.alliedhpi.com  
International offices in Europe, Asia, and South America

U.S.  
Toll Free: 800/444-3940  
Fax: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
Fax: 800/246-6201

International  
Telephone: 314/268-1683  
Fax: 314/771-5183